Pledget

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EN Instructions For Use (05)

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ES Modo de uso (14)

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SK Návod na použitie (26)

Pledget

Pledgets

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Pledgets

Description

Pledgets are pre-cut pieces from a sheet of either polytetrafluorethylene (PTFE) or expanded polytetrafluorethylene (ePTFE).

Indications

Pledgets are used as sterile, implantable devices for general reinforcement of non-absorbable surgical suture.

Contraindications

None known.

Warnings

Prior to using pledgets users should be familiar with surgical procedures and techniques involving pledgets and select the appropriate pledget-version according to the surgical requirements. Depending on the individual case, a non-appropriate fixation may lead to either leakages or necrosis and may support granuloma formation and/or pledget migration.

Do not expose PTFE products to temperatures above 260°C. PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.

Do not combine with absorbable suture.

This product is designed for single use only and may not be reused due to hygienic hazards.

Precautions

The handling of pledgets requires utmost care to avoid damage of the material.

When using forceps or needle holders, the material can easily be damaged through strong pushing or crushing.

When performing an autonomous trimming of the pledgets, it should be ensured that this results in an accurate and clean edge.

Dispose pledgets which were simply removed from the set, even when not actually used.

In case of an explantation pledgets may be disposed in accordance with common hygiene criteria.

Adverse reactions

On rare occasions pledget migration may occur which - in the worst case - can lead to vascular occlusion.

Like all foreign bodies, pledgets can reinforce an existing infection. A prolonged exposure to saline solutions such as those found in the urine or gall-bladder tract may result in calcification.

The product may also cause granuloma formation, dehiscence of wounds, wound infection as well as temporary local tissue inflammation.

Especially when using soft felt, filament fibers may detach occasionally.

Particularly when using larger pledgets, the scarring of the surrounding tissue may lead to shrinkage of the pledget.

How supplied

Pledgets are available in in several material qualities, sizes and shapes to suit the requirements of different procedures.

Symbols

Article number

Batch number

Use until (YYYY-MM)

Do not reuse

Do not re-sterilize

See instructions for use

Store dry

Do not expose to sunlight

Do not use if package is damaged

CE-mark and identification number
of notified body

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Santec Medicalprodukte GmbH