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## Optilene®

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Optilene® en Longpack to minimize the memory effect of the suture.

Optilene® with Premipatch® PTFE pledgets made from 100% polytetrafluoroethylene (PTFE)

- Optilene® size USP 5/0 to USP 8/0 sutures for cardiovascular indications with almost square bodied CV Pass needles with a thicker layer of silicone coating.

For a more detailed description of the product range see the respective product catalogue.

### Indications for use

Optilene® is indicated for soft tissue approximation and/or ligation in general surgery, when surgical practice requires the use of a non-

absorbable suture material. Optilene® sutures are also for use in cardiovascular surgery and neurosurgical procedures, where the low thrombogenic properties of polypropylene have a positive influence on wound healing. Other indications for Optilene® include skin closure, microsurgery and ophthalmology.

The Premipatch® PTFE pledgets are intended for the use as pads between the suture and the tissue surface in order to distribute the load of the suture over a bigger tissue area.

CV Pass needles are indicated for cardiovascular surgery.

**Contraindications**
The implantation of this suture is contraindicated in patients with known hypersensitivities or allergies towards its components (polypropylene, polyethylene, copper phthalocyanine, stainless steel and silicone).

Use of cutting needles is contraindicated for suturing parenchymal organs or other friable tissues where a cutting needle is not indicated.

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Upon the employment of Optilene®, the suture elicits a mild acute inflammatory reaction in the tissue followed by an encapsulation of the suture by fibrous connective tissue. Optilene® sutures are not absorbed by the human body and are not known to lose a significant amount of its tensile strength after implantation.

**Warnings**
Optilene® and all its special presentations must not be re-sterilized. Open unused or damaged packs shall be discarded. Do not reuse: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

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Title: FOLLETO OPTILENE V.5 Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)  
Title: SUT-PR-ES04-Técnico Sistemas Impresión  
Date: Monday, 17 January 2022, 11:54 W. Europe Daylight Time  
Meaning: Document signed as Author

---

---

UserName: García, Guillermo (garcgues)  
Title: SU-PR-ES04-Technical support manager  
Date: Monday, 17 January 2022, 14:23 W. Europe Daylight Time  
Meaning: Precheck of Document

---

---

UserName: Piera, Maria Eva (piermaes)  
Title: Senior Regulatory Affairs  
Date: Tuesday, 18 January 2022, 09:06 W. Europe Daylight Time  
Meaning: Approve Document

---

---